IMPLANT LINE

IFU_IL0008_ML_REV.C_2015-11-06

Abformpfosten

2

DE

Implant pick up

3

EN

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4

FR

Transfer per impronta

5

IT

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6

ES

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7

PT

Hersteller
Manufacturer
Fabricant
Produttore
Fabricante
Fabricante

Medentika® GmbH
Hammweg 8-10
DE-76549 Hügelsheim/Germany
Tel: +49-7229-69912-0
Fax: +49-7229-69912-20
E-Mail: info@medentika.de
www.medentika.de

Microcone® / Quattrocone® Abformpfosten

1. Indikation:

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C €0483

1. Indication:

The implant pick up is used around the implant and its exact position by taking an impression to be transferred to the model.

2. Contraindications:

none known

3. Warnings:

The instructions for use must be read before using the prosthetic and auxiliary components. These components must only be used in accordance with the indication under general accepted rules for surgical and dental treatment and as specified by workplace safety and accident prevention regulations. The components are intended for medical and dental use in combination with the Microcone® / Quattrocone® Implant System only. They must only be used in perfect condition. If the indication or the nature of use is not clear, do not use them until all points have been clarified. For further information, please refer to our Surgery Manual. We warrant that our products are of good quality under our terms and conditions of sale and supply. The following descriptions are not sufficient to ensure that dentists and surgeons who are not experienced with implant-prosthetic treatment can use the components correctly. We therefore recommend instruction in use by an experienced user. The user must be familiar with dental surgery and prosthetics, including diagnostics and preoperative planning and laboratory procedures. Because we have no control over the use of the product, no liability for any damages caused by its use can be accepted. The dentist or surgeon takes sole responsibility. The following complications have been observed occasionally:

- Aspiration or swallowing of components used in the mouth of the patient.
- Fracture of the straining screw in the abutment as a result of excessive torque.

4. Precautions:

The following precautions are to be met prior to or during treatment.

- Prior to each procedure it must be ensured that all necessary components, instruments and materials are available in the requested quantities.
- Always wear protective clothing for your own safety.
- Position the patient so that the danger of aspiration of components is minimized. All components that are used intra-orally must be secured to prevent aspiration or swallowing.

5. Adverse Reactions:

none known

6. Application:

Implant type	Implant connection	Compatible components
Microcone	NI (Narrow Interface) D 3,0 RI (Regular Interface) D 3,5-5,0	Identification: NI Identification: RI
Quattrocone	RI (Regular Interface) D 3,5/4,3	Identification: RI
Quattrocone30	AI (Angulated Interface) D 4,3	Identification: AI

Impression taking

Open tray technique/Pick up technique

Remove gingiva formers, spray tapered connection with air/water spray and dry them. Select implant pick up, position in one of the four possible index positions, check correct seating and fix manually with transfer screw. The hex is solely for loosening the screw. Impression taking with open tray. After hardening loosen transfer screw and remove the impression with the implant pick up.

Closed tray technique/repositioning technique

CAUTION: Repositioning technique is contraindicated where the implants have axial divergences.

Remove gingiva former, spray tapered connection with air/water spray and dry them. Place the closed tray implant pick-up in one of the four indexing positions, check that it is seated correctly and tighten with the SW 1.26 mm screwdriver. Then put the repositioning cap to the end position on the implant pick up. The impression is made with soft silicone, hydrocolloid, or polyethylene compounds. After hardening remove the impression and remove implant pick up from the mouth. Mount implant pick up on the Lab Analog and then reposition the implant pick up in the impression. Make sure that the impression material is not pinched.

7. Cleaning/Sterilization:

The implant pick up is delivered in a non-sterile state in the appropriate packaging. The products has to be cleaned and sterilized according to following details:

Cleaning:

Automatic cleaning procedure based on Vario TD program:

- Cold water rinse 1 min.
- Cleaning at 55°C/131°F (± 2°C/3.6°F) for 5 min. with cleaning agent neodisher® MediClean (0.5% v/v)
- Neutralization with cold water for 2 min.
- Cold water rinse 1 min.

Sterilization:

Seal in a suitable autoclave foil.

Steam sterilization with 3x fractionated pre-vacuum: 132°C (270°F) for 4 minutes, drying time: 20 minutes

Signs and Symbols:

- For single use only
- LOT number
- Reference number
- Manufacturer
- Attend instruction for use
- Non sterile

Products indicated with ® are registered brand names of the manufacturers.

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