




PreFace abutments

Straight

Titanium Grade 5 ELI



Indications

- Single tooth or multiple unit implant restorations
- Cement-retained restorations

1. Impression (Dentist)

- Remove the healing abutment or temporary restoration from the implant
- Spray implant interface with air/water spray and dry them
- Select Implant Pick-up and place it into one of the possible index positions of the implant
- Check correct seating and fix it manually
- Take an implant level impression with suitable impression material
- Opposing jaw information and shade are included with the impression and Implant Pick-Up when sent to the dental laboratory



2. Laboratory procedures (Dental technician)

- Make sure that the impression material is free of cracks and bubbles
- Mount corresponding Lab Analog on the Implant Pick-Up and prepare a master hard-plaster model with a gingival mask that can be scanned and removed
- Insert the scanbody into the Lab Analog and screw it in with the corresponding screw.
- Select the appropriate implant and PreFace abutment in the Medentika library and scan the model with inserted scanbody.
- Scan the gingival mask
- Construct the custom-made emergence profile and the oral part of the abutment in accordance with the anatomical specifications.

Ensure that the minimum wall thickness of the PreFace abutment is not less than 0,5mm and the angle does not exceed 30°. The diameter and height of the emergence profile should not exceed 13mm and 6mm, respectively. The total height of the abutment should be between 4,0 and 15,0 mm.

- Send the data file to the Straumann milling center.

3. Milling of the Preface abutment (Straumann milling center)

- Upload the data file and check construction recommendations regarding wall thickness, angulation, emergence profile and length.
- Select the correct PreFace abutment and mount it into the PreFace abutment holder
- Mill the PreFace abutment according to the transferred data file.
- Send the customized Preface abutment to the dental technician



4. Abutment checking (Dental technician)

Note: *For strength and fit reasons, never modify the implant interface of the PreFace abutment*

- Insert customized Preface abutment into the Lab Analog of the hard-plaster model
- Check correct construction, fit to the Lab Analog and anatomical situation.

5. Abutment connection (Dentist)

Note: *For strength and fit reasons, never modify the implant interface of the PreFace abutment*

- Ensure that the implant platform is free from any soft-tissue or bone remnants.
- Position the abutment into the implant and secure the screw in the implant using the suitable screwdriver.
- A radiograph can be helpful to confirm accurate seating of the abutment
- Tighten the abutment screw to the torque that is recommended in the instruction for use of the abutment using a torque ratchet.

6. Cementation of final restoration (Dentist)

- Gently seat the restoration on the abutment and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Fill the screw access channel with a block-out material to preserve abutment screw access.
- Cement the restoration using permanent cement according to the instructions for use of the cement manufacturers