## ABUTMENT LINE

IFU\_AL0012\_ML\_REV.C\_2015-02-04

|                                   |    |    |
|-----------------------------------|----|----|
| Drehmomentratsche                 | 2  | DE |
| Torque control ratchet            | 4  | EN |
| Clé à cliquet dynamométrique      | 5  | FR |
| Chiave dinamometrica a cricchetto | 7  | IT |
| Carraca dinamométrica             | 9  | ES |
| Chave de torque                   | 11 | PT |

Hersteller  
Manufacturer  
Fabricant  
Produttore  
Fabricante  
Fabricante

Medentika® GmbH  
Hammweg 8-10  
DE-76549 Hügelsheim/Germany  
Tel: +49-7229-69912-0  
Fax: +49-7229-69912-20  
E-Mail: [info@medentika.de](mailto:info@medentika.de)  
[www.medentika.de](http://www.medentika.de)

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## Instructions for use torque control ratchet M 12

# CE0483

**1. Use of the torque control ratchet**

The torque ratchet consists of the ratchet head (Fig. 1-a), ratchet handle (Fig. 1-b) and adjusting screw (Fig. 1-c). Medentika placement instruments manual/ratchet (compatible with the Straumann system) can be used in the ratchet head (Fig. 1-a).


Fig. 1

The ratchet head is lockable to block the torque control (see Section 1.1). When the ratchet head is not locked the torque can be adjusted progressively from 10 to 40 Ncm (see Section 1.2).

Before initial use and after each use on a patient, the ratchet must be dismantled, cleaned, dried and **sterilised immediately before further use** (see Section 3).


Fig. 2

Fig. 3

Fig. 4

Fig. 5

**1.1 Locking the ratchet**

The ratchet head is supplied unlocked. If the ratchet head should be locked, turn the adjusting screw to lower than 10 Ncm (Fig. 4), so that the spring is fully released. Then pull the ratchet handle out of the ratchet head (Fig. 2), turn the ratchet handle clockwise by 90° and lock it in place in the ratchet head (Fig. 3). Then set the torque to 20 Ncm to ensure that the ratchet remains locked.

**Attention: Adhere to the max. permitted loading of 100 Ncm in the rigid condition!**

**1.2 Setting the torque function**

First, ensure that the ratchet is not locked. Turn the adjusting screw until the required torque is reached. The torque to be set is reached when the respective marking ring of the adjusting screw on the end of the handle is just visible (Fig. 4).

When loaded, the ratchet head now bends once the set torque is reached (Fig. 5). It is essential not to turn any further, as otherwise there is the risk of greatly exceeding the torque.

**Important note: When the ratchet is not in use the adjusting screw should be turned back until the spring is fully released.**

**1.3 Changing instruments** (see at 2.1.1 and 2.2.5)

**2. Maintenance of the torque ratchet**
**2.1 Dismantling**

- 2.1.1 Remove the manual/ratchet insertion tool (Fig. 6-e) from the ratchet head (Fig. 6-a).
- 2.1.2 Fully unscrew the adjusting screw (Fig. 6-d) from the ratchet handle (Fig. 6-b).
- 2.1.3 Pull the ratchet handle (Fig. 6-b) out of the ratchet head (Fig. 6-a) and remove the spring (Fig. 6-c).
- 2.1.4 Disinfection, cleaning and drying the elements (see Section 3). The ratchet should then be assembled before sterilisation.


Fig. 6

**2.2 Assembly**

- 2.2.1 All elements should be disinfected, cleaned and completely dry before assembly. Lubricate (see Fig. 6) the areas marked with ↗. This reduces the metal to metal friction and prevents friction corrosion. Only instrument oils (e.g. paraffin-based oil spray) with approval for steam sterilisation and tested biocompatibility should be used.
- 2.2.2 Place the ratchet handle (Fig. 6-b) over the threaded rod on the ratchet head (Fig. 6-a).
- 2.2.3 Slide the spring (Fig. 6-c) into the ratchet handle (Fig. 6-b).
- 2.2.4 Wind the adjusting screw (Fig. 6-d) into the ratchet handle (Fig. 6-b) until the resistance of the spring is felt.
- 2.2.5 Insert the required insertion tool manual/ratchet (Fig. 6-e) in the ratchet head (Fig. 6-a).

**Important note: Before sterilisation the ratchet should be checked for correct functioning in both directions.**

**3. General preparation**
**3.1 Preliminary conditioning**

Place the ratchet in a cleaning/disinfection agent immediately after use. This is for personal safety and prevents contamination drying onto the ratchet. The cleaning agent manufacturer's instructions regarding concentration, immersion time and temperature should be strictly observed to avoid corrosion and discoloration.

The cleaning and disinfection agent used during conditioning is only for your personal protection and – after cleaning – cannot replace the disinfection stage to be completed later.

**3.2 Cleaning and disinfection**

The ratchet must be dismantled for cleaning and disinfection (see at 2.1).

Cleaning and disinfection can be performed automatically with the aid of a thermodisinfector (adhere to the manufacturer's instructions) or manually using suitable cleaning and disinfection solutions.

Use of an enzymatic cleaning agent with an almost neutral pH value is recommended. Do not use cleaning and disinfection agents that contain halogens (chlorine, iodine, bromine, fluorine), oxidation agents (H<sub>2</sub>O<sub>2</sub>), organic solvents, heavy metal salts, aldehydes or aromatic/halogenated hydrocarbons.

Accretion, contamination etc. should be thoroughly removed using a soft brush (never use metal brushes, abrasive cleaners, steel wool or apply excessive force).

**Note:** The manufacturer's instructions should be observed when using an ultrasonic cleaner for disinfection and cleaning.

The following stages are recommended when performing the **automatic cleaning/disinfection process in a thermodisinfector**:

**Cleaning**

- Rinse for 1 minute with cold water
- Clean for 5 minutes using a suitable cleaning agent at 55°C (± 2°C) (e.g. neodisher MediClean (0.5%))
- Neutralise for 2 minutes in cold water
- Rinse for 1 minute with cold water

**Disinfection**

- 5 minutes thermal disinfection at 93°C using demineralised water

The following stages are recommended when performing the **manual cleaning/disinfection process**:

**Cleaning**

- Mechanical cleaning of the entire surface, particularly of the hollow sections, gaps and drill holes under running water (drinking water, 30 ± 3°C) with the aid of a soft brush until all residue has been removed (20 seconds minimum per component).
- Clean for 5 minutes in an ultrasonic cleaner (minimum frequency: 35 kHz) using a suitable cleaning solution (e.g. 2% MediClean®, Dr. Weigert, in drinking water).
- Rinse the entire surface, particularly hollow sections, gaps and drill holes under running water (drinking water, 30 ± 3°C) for a minimum of 10 seconds per component. Moveable parts must be moved when rinsing.
- Rinse with demineralised water for a minimum of 5 seconds per component.

**Disinfection**

- Immerse the components for 30 minutes in disinfectant solution (e.g. 2% neodisher® SeptoMED, Dr. Weigert, in drinking water at 20 ± 2°C) and move the components in the solution to remove air bubbles on the surface.
- Immerse for 1 minute in cold, demineralised water.
- Thoroughly rinse the entire surface, particularly hollow sections, gaps and drill holes under running water (drinking water, 30 ± 3°C) to remove the disinfectant solution residue. Moveable parts must be moved when rinsing.
- Rinse with demineralised water for a minimum of 5 seconds per component.

**3.3 Sterilisation**

The ratchet should be sterilised assembled and unloaded, whereby only thoroughly cleaned and disinfected components without any signs of corrosion may be assembled and sterilised. It is recommended to complete the sterilisation using a vacuum autoclave. Pack the surgical tray with suitable sterilisation supports and sterilise using 3x fractional prevacuum for 4 minutes at 132°C (270°F). Drying time is 20 minutes.

Products marked with ® are registered trademarks of the respective manufacturer.

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